

DISPOSABLE FACE MASK

DESCRIPTION

- **Product Classification:**
 - Respiratory Protection - Disposable face mask
- **Construction:**
 - Triple-layer construction, ear loops
 - 100% non-woven polypropylene inner and outer mask, 100% melt-blown polypropylene interior material; Spandex/polyester blend ear loops, adjustable nose clip
- **Color(s):** Blue/white
- **Unit of Measure:** Each (50 each minimum order)
- **Packaging:** 50 each per box, 2000 each per case



FEATURES & BENEFITS

- Disposable face masks shield the mouth and nose from large-particle droplets, splashes, sprays, or splatters
- Triple-layer construction
- Non-woven material
- Comfortable stretch polyester ear loops
- Adjustable nose clip prevents eyewear from fogging
- Please refer to General Information for notes on proper usage

SERVICE LIFE

The service life of this product will vary depending on the application and type of care administered. This product should be replaced immediately once a critical or major defect (as defined below) is noted.

- **Critical Defects:** Holes/tears
- **Major Defects:** Dimensional or weight nonconformance, discoloration, other fabric defects

GENERAL INFORMATION

Notes on Proper Usage: If worn properly, a mask is meant to help block large-particle droplets, splashes, sprays, or splatter that may contain germs (viruses and bacteria), keeping it from reaching your mouth and nose. Masks may also help reduce exposure of your saliva and respiratory secretions to others.

While a mask may be effective in blocking splashes and large-particle droplets, a face mask, by design, does not filter or block very small particles in the air that may be transmitted by coughs, sneezes, or certain medical procedures. Masks also do not provide complete protection from germs and other contaminants because of the loose fit between the surface of the face mask and your face.

Masks are not intended to be used more than once. If your mask is damaged or soiled, or if breathing through the mask becomes difficult, you should remove the mask, discard it safely, and replace it with a new one. To safely discard your mask, place it in a plastic bag and put it in the trash. Wash your hands after handling the used mask.

These masks do not supply oxygen.

These masks are manufactured in a sterile environment, but are not certified for surgical settings.

Packaging Note: The size, material, and artwork on the packaging may vary. All MM0 masks will meet the above specifications.